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Manual for Cryogen Usage

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Cryogenics and Instrumental Analysis Division

Natural Science Center for Basic Research and Development

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A person who will transfer the liq. He must reserve the transfer time according to the table for cryogen usage. When you will transfer the liq. He, you should put the usage card in the box, and write necessary matters on the notebook.

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A beginner must transfer the cryogen with an experienced person. When any experienced persons are not in laboratory, you should receive a technician's guidance.

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User must recover the He gas. If user will do an experiment in which He gas can not be recovered, you should take steps for an extra charge before the experiment.

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Please refer to "Lending of Cryogen Vessels" (p.42) in this text. User should book the lending and write necessary matters on the note book in the control room (I-101). Please put the seal of a person responsible for expense burden.

### 4. 密閉型液体窒素容器の検査 Examination of a Closed Liq. N<sub>2</sub> Vessel

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The examination of a closed Liq. N<sub>2</sub> vessel (high-pressure vessel) in a certain period is required according to the High Pressure Gas Safety Act. Please refer to "Examination of a Closed liq. N<sub>2</sub> Vessel" (p.41) in this text, and then undergo the examination. User should book the lending and write necessary matters on the note book in the control room (I-101). Please put the seal of a person responsible for expense burden.

### 5. ヘリウム回収管の新設と利用 New Construction of a He Gas Recovery Line and its Usage

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Please refer to "New Construction of a He Gas Recovery Line and its Usage" (p.43, 44) in this text. Before new construction of a He gas recovery line, please contact us, and receive a technician's guidance. After the construction, please submit a report for the He gas-tight test. Please refer to "Note of the New Construction of a He Gas Recovery Line"(p.43) and "Leak Test after the New Construction of a He Gas Recovery Line and its Usage" (p.43) in this text.

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Liquid Helium (He) and liquid Nitrogen (N<sub>2</sub>) regarded as high-pressure gases for High Pressure Gas Safety Law in Japan are dangerous if they are used incorrectly. Furthermore, Helium is a limited natural resource imported from USA mainly. Therefore, we should understand these properties of liquid He and N<sub>2</sub>, and use them efficiently and safely.

At first, fundamental points in the usage of cryogenes are mentioned in the section 2. In sections 3 and 4, the way of transfer of cryogenes is explained in more details, and some manuals for usage of cryogenes are listed up in section 4. But, because these are minimal know-how, persons carried out low-temperature experiments should learn some detailed textbooks for low-temperature techniques.

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Unit of pressure is "Pa" in International System of Units. There are many units of pressure such as kgf/cm<sup>2</sup>, psi, Torr, which are used conventionally. Especially, "psi" is used in the pressure gauge of liquid He vessel. In High Pressure Gas Safety Act, gauge pressure is used. The gauge pressure is based on the atmospheric pressure (=zero), usually used in an engineering field. In physical and chemical fields, absolute pressure, which is based on vacuum, is used.

Pa(=N/m <sup>2</sup> )	Torr(=mmHg)	kg/cm <sup>2</sup>	psi(lb/in <sup>2</sup> )	atm(気圧)	水柱 (m)
1	7.50062x10 <sup>-3</sup>	1.020x10 <sup>-5</sup>	1.4504x10 <sup>-4</sup>	9.869x10 <sup>-6</sup>	1.0206x10 <sup>-4</sup>
133.322	1	1.360x10 <sup>-3</sup>	1.9337x10 <sup>-2</sup>	1.316x10 <sup>-3</sup>	1.3607x10 <sup>-2</sup>
9.80665x10 <sup>4</sup>	735.559	1	14.2233	0.967841	10.0090
6.89476x10 <sup>3</sup>	51.7149	7.031x10 <sup>-2</sup>	1	6.805x10 <sup>-2</sup>	0.703704
1.01325x10 <sup>5</sup>	760	1.03323	14.6959	1	10.3416
9.79781x10 <sup>3</sup>	73.4896	9.991x10 <sup>-2</sup>	1.42105	9.670x10 <sup>-2</sup>	1

1 bar=10<sup>5</sup> Pa =1.0197 kg/cm<sup>2</sup>=0.98692atm=14.504 psi, psi: pound per square inch  
lb: libra(重量ポンド), 1 (米国 ton)/ft<sup>2</sup>= 0.94508 atm, 1 Pa = 1 N/m<sup>2</sup>=10 dyn/cm<sup>2</sup>

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The High Pressure Gas Safety Act was enacted to prevent accidents and disasters caused by high pressure gases. As mentioned below, high pressure gases become very dangerous if we make a wrong use of those for the production, storage, and consumption. Therefore, persons who produce the high pressure gases as well as who use those (students and staff of university) must comply with laws and other related acts, and pay attention to preventing accidents and disasters and securing public safety.

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Among high pressure gases, compressed gas is currently gas, and liquefied gas is currently liquid. High pressure gases defined by High Pressure Safety Act is not only specified materials but also compressed gas having specified states and liquefied gas having specified properties. (High Pressure Safety Act, Article 2).

1. Compressed gas, the pressure of which is not less than 1 megapascal at a temperature of 35 degrees Celsius (except compressed acetylene gas in both cases);
2. Compressed acetylene gas, the pressure of which is not less than 0.2 megapascal at a temperature of 15 degrees Celsius;
3. Liquefied gas, liquefied gas, the temperature of which is 35 degrees Celsius or less in the case that the pressure is 0.2 megapascal;
4. In addition to what is listed in the preceding item, those liquefied gases, the pressure of which exceeds zero Pascal at a temperature of 35 degrees Celsius, and which, inclusive of liquefied hydrogen cyanide and liquefied methyl-bromide, are specified by a Cabinet Order.

The pressure of the liquefied gas is identified by the saturated vapor pressure. The temperatures for L. N<sub>2</sub> and L. He at which these saturated vapor pressures is more than 0.2 MPa are about -

189 °C and -268 °C, respectively. Therefore, both L. N<sub>2</sub> and L. He are high pressure gases.

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#### Production, Storage and Consumption of High Pressure Gases

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The production of the high pressure gas is to make artificially the state of a high pressure gas be means of compression, liquefaction and some other means. We need an application and notification to the competent authorities according to the amount of a processing gas.

e.g.: "Production of the high pressure" is making a high pressure gas into the higher pressure one. To the contrary, it is making the high pressure gas into the lower pressure one. (The case that the pressure of an outlet becomes more than 0.2 MPa for the transfer of the liquefied gas, and so on.)

The storage of the high pressure gas is keeping the state of the high pressure gas of more than a certain amount at a fixed place. We need an application and notification to the competent authorities according to the storing amount of the gas.

e.g.: Cold evaporator (CE) for a L. N<sub>2</sub>

The consumption of the high pressure gas is making the high pressure gas into the state of not being the high pressure gases by means of a relief valve and so on.

e.g.: Transfer of a high pressure gas in a container to an open type container. Putting the gas of a bomb into a balloon through a relief valve (regulator).

There is no regulation for the consumption of a small amount of gases, but regulation for a technical method for the consumption. We need a notification to the competent authorities in advance.

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#### Cryogenics Section and Safety Training

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
The history of the Cryogenic section of Cryogenics and Instrumental Analysis Division in Natural Science Center for Basic Research and Development (N-BARD) started since the establishment of the Low-temperature room in Faculty of Science in 1966. The Cryogenic section has existed up to now through Cryogenic Center established in 1988. Now, the staffs of the Cryogenic section produce a liq. He and store liq. N<sub>2</sub>, provide shared use apparatuses and experimental rooms, and enlightening activities of low temperature science as Class 1 Producer specified by paragraph (1) of Article 5 of High Pressure Gas Safety Act.

Class 1 Producer has been required the subjecting our Designated Facilities, periodically (once a year), to a safety inspection conducted by the prefectural governor, and drawing up a Hazard Prevention Rule and submitting notification report thereof to the prefectural governor, preparing a Safety Training Program and implementation of the Training, conducting self inspections periodically for safety and daily inspections. At Cryogenic section, the staffs make efforts to supply Liq. He and Liq. N<sub>2</sub> stably by means of the maintenance of manufacturing and storage facilities and voluntary operational safety activities etc.

At Cryogenic section, the staffs conduct the work for the safety of the production of cryogenes by organizing the safety management texture as below. The Safety Training for workers is regulated in Article 29 of the Hazard Prevention Rule, and Article 27 of High Pressure Gas Safety Act. The cryogen usage techniques and safety lecture doing by cryogenic section is carried out on the basis of their Articles as a part of the Safety Training.

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Safety management texture in Cryogenics section

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Table 1. shows the summary of the physical properties of some typical cryogens. These possess low boiling point, relatively lower latent heat, large volume ratio of the gas and liquid. In using the cryogen, there exist possibilities of **accidents of frostbite, lack of oxygen (suffocation), explosion, and others**. It is important to understand their properties very well and handle the cryogens.

Table 1. 寒剤の性質 Cryogen properties

液体 Liquid	分子量 Molecular weight	沸点 (K) Boiling point	密度 (g/cm <sup>3</sup> ) <sup>†</sup> Density	潜熱 (cal/cm <sup>3</sup> ) Latent heat	容積比* Volume ratio
O <sub>2</sub>	32.0	90.19	1.144	58.2	800
N <sub>2</sub>	28.01	77.35	0.808	38.5	647
<sup>4</sup> He	4.003	4.216	0.125	0.62	700
<sup>3</sup> He	3.000	3.1905	0.059	0.12	440

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Volume ratio: the ratio between the volume of the liquid at the boiling point and that of the gas at 0°C and 1 atm.

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In the handling the cryogen, your hands and foots may be frostbitten by causes mentioned below.

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When the bare hands or body are exposed to a cryogen or cold gas.

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When you touch a cold metal such as the surface on a delivery tube. When a small amount of liq. N<sub>2</sub> drops on your hands, your hands protected by the N<sub>2</sub> gas is not frostbitten. However, when a cold metal cooled by the liq. N<sub>2</sub> is touched by bare hands, the hands and metal sticks together because of the no protection of the gas, and then, the hands are heavily frostbitten. In the cold metal cooled by the liq. He, there becomes a serious situation. (See reference materials p. 55)

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When you touch cryogenics directly by your wet hands or feet.

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You should use leather gloves or thermal insulating gloves in the handling the cryogen. Do not use cotton or woolen work gloves in which liquid N<sub>2</sub> sinks.

Fig. 2. 断熱手袋, 革製手袋と軍手 Thermal insulating gloves, leather gloves and cotton work gloves

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Please brush away the Liq. N<sub>2</sub> quickly when the liq. N<sub>2</sub> splashes.

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You should use goggles to prevent the liq. N<sub>2</sub> from getting into the eye when the Liq. N<sub>2</sub> splashes.

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You should soak the frostbitten skin in a hot water. Do not use a drier. After that, go to a hospital.

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When your hand stick to a cold metal, you warm up the metal up to the temperature of tepid water, and then, take your hand off the metal slowly. After that, you should go to a hospital.

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When the Liq. N<sub>2</sub> enters your eyes, you should wash your eyes by clean water.

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Vaporization of liquid N<sub>2</sub> will increase its volume and decrease the concentration of O<sub>2</sub> in the room. **When the cryogenics are used in the closed room, we must ventilate the room.**

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The rapid oxygen-shortage is caused by breathing a gas without oxygen. **When you come in the room in which the oxygen is 0%, you fall down in one breathing, and have no chance of resuscitation.** To breathe near the evaporating gas of liq. N<sub>2</sub> is a suicidal act.

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The gradual oxygen-shortage is caused slowly by breathing an air in which the oxygen concentration decreases gradually. (See Table 2.)

Table 2. 酸素濃度低下(緩急性酸欠)による人体への影響  
Effects on the Human Body due to Gradual Oxygen-shortage

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[酸欠事故の被害拡大の要因] Factors of the Accident in Oxygen-shortage

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Secondary Disaster: There is a high possibility that a rescuer coming into the room without oxygen suffer damage because he can not find why a person falls down in the room without smoke or irritating smell.

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Lack of Judgment Ability: A person coming into the room without oxygen can not understand that his headache, nausea, and coma are attributed to the lack of oxygen.

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Sharp Drop in Resuscitation Opportunity: A chance of resuscitation decreases rapidly in 3~5 minutes after the stop of breathing, and becomes zero above 7 minutes (Fig. 3). Even if a person is resuscitated, he has sometimes after effects of disease such as speech impediment, motility disturbance, narrowing of visual field, hallucination, neurosis.

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A rapid treatment preventing the secondary disaster is required.

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When you use liq. N<sub>2</sub> in a closed room, you should release the N<sub>2</sub> gas out the room by using a ventilating fan or exhaust pipe. (See reference materials p.54)

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When you take a vessel of liq. N<sub>2</sub> by car, you must open the windows.

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When liq. N<sub>2</sub> or He is transported by elevator, the vessels should travel alone in the elevator, without any people accompanying them inside the elevator. (See reference materials pp.45, 46)

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When the inside of a room seems to be in the lack of oxygen, you should ventilate the room with a fresh air, and you may enter the room after you find that the oxygen concentration becomes 20% or more. You can examine the oxygen concentration by an oxygen tester borrowed from Cryogenics section.

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Liq. He and liq. N<sub>2</sub> boil even in the vessels because the heat enters from outside. When the liq. He and N<sub>2</sub> are transformed into gases at 273 K and 1 atm, their volumes expand up to 700 and 647 times, respectively. When cryogen vessels are sealed, there is a danger of explosion because the pressure in the vessels becomes very high. **Do not seal cryogen vessels. (See Table 1.)**

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The vent of the vessel containing a cryogen is closed by ice.

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All of valves in the sealed vessel are closed.

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You should mount check valve in the vent by preventing air entering into the vessels.

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When you use relatively large sealed vessels with some relief valves, please read those manuals and understand how to use them very well.

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### —The Neck in the Vessel is Weak Mechanically

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The neck in the vessel is weak mechanically. Therefore, do not knock the vessel down, and put the large vibration.

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### —Notes for Using the Vessels and Transfer of the Cryogenics

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If there are some contaminations and water in the vessel, we should remove them before the transfer of cryogenics.

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It is not good that cryogen is poured suddenly into a vessel at room temperature.

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There is a possibility that the large vessel over 100 L stumbles over a stone.

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Please open the cap of the closed vessel after the inner pressure is released.

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Do not use a safety valve as a release valve.

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**Do not leave the place where we transfer cryogenics.** In the transfer of cryogenics, we have possibility to be in unexpected accidents. We always need to watch circumstances carefully.

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**Do not handle anything except transfer valve in liquid N<sub>2</sub> tank.** There are some relief and safety valves in liquid N<sub>2</sub> tank. If a user handles these valves, there is possibility to increase the pressure

in the tank. **Please handle only transfer valve in liquid N<sub>2</sub> tank.** If you find some anomalies, please contact technicians in cryogenics division.

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**Close the valve tightly after the transfer of cryogenics.** The valve is frozen because the cryogenics come out little by little.

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### —Pay Attention to Safety and Suppression of Cryogen Losses in the Transfer of Cryogenics

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In the transfer of cryogenics, we should carry out precooling by the cold gas. At first, we put the transfer tube into the bottom of the vessel, and cool slowly the tube and inner part of the vessel.

寒 剤 Cryogenics	He		N <sub>2</sub>
冷却温度 Cool down temperature (K)	4.2		77.3
初めの温度 Initial temperature (K)	300	77	300
潜熱だけを利用 Case for using the latent heat of the gas	31	2.2	0.46
潜熱とガスエンタルピーを利用 Case for using the latent heat and enthalpy of the gas	0.79	0.15	0.29

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When liquid N<sub>2</sub> is exposed to the air, the oxygen is condensed in the liquid N<sub>2</sub>. We have the hazard of explosion when flammable gases or organic compounds come into the liquid N<sub>2</sub> with the oxygen. In order to avoid a back current of the air, please use Bunsen valve or bubbler.

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### —Do not Handle Roughly the Transfer Tube for Liq. He.

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As shown in Fig. 10, the transfer tube for Liq. He has a double tube structure, which has an insulation vacuum. If it is handled roughly, it can not be used because the spacer is broken and the inner and outer tubes are touched each other or the vacuum gets worse. The outer diameter of the tube in the center is 12mm.

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### —Do not Put an Air into the Recovery Line for He Gas

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The connections between the He recovery line and vessel or cryostat should be made tight. If the vessel would be disconnected from the recovery line, we have the hazard of high pressure in the vessel because an ice blocks the opening.

When you pump the He gas in the recovery line, the cryostat and pump should be made airtight completely.

### 3. 寒剤汲み出しの手順と諸注意 —Notes and Procedures for the Transfer of Cryogens

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##### —Data for cryogens usage, Booking procedure, service area—

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Liq. N<sub>2</sub> is not supplied between 11:30 and 13:00. The transfer of Liq. N<sub>2</sub> should be finished until 16:30.

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For Liq. He, a reservation alternation and added reservation after deadline can not be allowed.

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When you will transfer the liq. He and N<sub>2</sub>, you should put the usage card in the box, and write necessary matters on the notebook.

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#### Notes for filling in the cryogen usage card of Liq. N<sub>2</sub> and He

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The cryogen usage card is a basic record for using a cryogen. It is very important document to calculate the user fees for cryogenes. The card is supplied in this center. Please note these points mentioned below for filling out.

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Procedures for the transfer of L. He are mentioned in the next section. Please follow the technician's guidance for a detailed procedures of the L.He transfer. Please refer to the section 3.1 for booking cryogen usage.

Procedures for the transfer and usage of L. He is listed below.

1. Pre-cooling of the vessel
2. Transfer of L. He to the vessel
3. Usage of L. He and recovery of the He gas

The operation 1. Is not needed if the L. He remains in the vessel, or the temperature in the vessel is down to L. N<sub>2</sub> temperature (77K). Please refer to the section 2.8, and 2.9 (p. 15) in the chapter 2 for the notes of L. He usage.

#### 3.3.1 大口の液体ヘリウム容器の予冷 Pre-cooling of the large volume vessel

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Recent large volume vessel is a gas shielded type, which is shielded thermally by the cool He gas evaporated in the vessel. For the pre-cooling of the vessel, you should put L. N<sub>2</sub> into the vessel before 3 days of the transfer of the L. He. Just before the transfer, you expel the L. N<sub>2</sub>, and replace the N<sub>2</sub> gas with He gas. To avoid a back current of the air into the vessel, please note the seal. After the L. He is transferred, it is preferable that the small amount of L. He remains in the vessel, because of the reduction of the cryogen losses in the transfer.

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Please do pre-cooling by next procedures when the temperature in the vessel is higher than L.N<sub>2</sub> temperature.

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How to put L. N<sub>2</sub> into the vessel (before 3 days of the transfer of the L. He)

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Please remove the quick-coupling on the top of the vessel and see visually that there are some contaminations and/or water in the bottom into the vessel or not. If there are something, please remove them.

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Put L. N<sub>2</sub> into the vessel from D. The injection amount is 1/3 of the vessel volume.

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Please confirm that the temperature in the vessel is L. N<sub>2</sub> temp.

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Put the SUS tube of 12 mm in outer diameter into the vessel from D to the bottom of the inside of the vessel.

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Remove the L. N<sub>2</sub> by putting the He gas (0.049 MPa) into the vessel from E. After removing the L. N<sub>2</sub>, remove SUS tube and close valve B.

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Confirm removing L. N<sub>2</sub> completely by seeing visually into the vessel.

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Put He gas of 0.02 MPa into the vessel from E.

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(Notice) Transfer of L. He should be canceled if L. N<sub>2</sub> remains in the vessel.

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#### How to transfer L. He to the vessel with a large volume

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Put the cryogen usage card to the box labeled "Liq. He".

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Write necessary matters such as date, time, department, name of a person responsible for expense burden, user name, volume of vessel, liquid He volume, and number of vessel on the notebook.

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Put the vessel on the table lifter and connect the flexible tube for recovery of G. He to the vessel. Confirm the open of the valve M and then open valve A. After that, remove the nut on the insertion hole and see the O-ring and O-ring presser.

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Take the lift down to the position in which the top of the transfer tube can be put into the insertion hole safely.

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Set the top of the transfer tube on the insertion hole. At that time, if there are some drops of water on the surface of the tube, wipe away them. Unclamp the nut on the insertion hole and open valve B and then put the transfer tube into the insertion hole. Take the lift up to the top as you hold the tube by hand not to strain and move the tube.

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How to transfer the L. He for large vessel.

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
Fig. 14. 液量測定  
Measurement of L. He volume


Fig. 15. 液量測定 (拡大図)  
Measurement of L. He volume (expansion)

Liquid level gauge is using Taconis oscillation which is a gas oscillation in a thin tube with a large temperature gradient between top and bottom. The change of frequency when an end of the thin tube closed another end contacts to a liquid surface allows you to find a height of L. He.

Way to determine the liquid volume follows next procedures. See that valve A opens. Remove the upper lid of Fig. 13 and open valve B. After that, put the liquid level gauge into the bottom of the vessel, in which the position is basis ( $h_0$ ) of the height of L. He. When you lift slowly the tube from that position, the sound of the vibrating plate on the top of the gauge becomes bigger at a position. Please determine the position ( $h_1$ ) in which the sound becomes bigger when the tube is lifted up. Measure the length of  $h_1-h_0$  as shown in Fig. 15. After that, remove the gauge from the vessel and close valve B. Put the upper lid of Fig. 13 and gauge back to the original place. Calculate the L. He volume using a conversion table with the vessel.

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#### Manual for usage of self-pressurized L. N<sub>2</sub> vessel (closed type vessel)

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There is self-pressurized vessel in some L. N<sub>2</sub> vessels with rather large volume. L. N<sub>2</sub> comes out from the vessel by increasing the pressure in the vessel to evaporate L. N<sub>2</sub> introducing to a rising pressure tube inside the vessel. This vessel is often used because other pressure source is not necessary. This vessel is a closed type and is regulated by High Pressure Gas Safety Act. Please be careful when handling the vessel. See Examination of a closed type Liq. N<sub>2</sub> Vessel (p.41).

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Close the valve A and B, open the valve C.

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After the pressure in the vessel is confirmed to be zero, remove the plug D, and then Put L. N<sub>2</sub> into the vessel.

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If valve C opens, close it. And then, open valve A.

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Open valve B little by little, and then G. N<sub>2</sub> comes out little by little. Hold the pressure of 0.2-0.5 kg/cm<sup>2</sup>≒0.02-0.049 MPa to control valve B. When you take out enough volume of L. N<sub>2</sub>, close valve A and B.

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(Notice)

- ✓ Use the vessel within the maximum filling pressure.
- ✓ Prevent that water drops attach the safety equipment such as relief valve and rupture disk and the equipment is frozen.

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Article 1 The purpose of this Act is to regulate the production, storage, sale, transportation and other matters related to the handling of high pressure gases, their consumption as well as the manufacture and handling of their containers and to encourage voluntary activities by private businesses and the High Pressure Gas Safety Institute of Japan for the safety of high pressure gases with the aim of securing public safety by preventing accidents and disasters caused by high pressure gases.

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(i) Compressed gas, the pressure (meaning gauge; the same shall apply hereinafter) of which is not less than 1 megapascal at its normal operating temperature and which is currently not less than 1 megapascal, or compressed gas, the pressure of which is not less than 1 megapascal at a temperature of 35 degrees Celsius (except compressed acetylene gas in both cases);

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(ii) Compressed acetylene gas, the pressure of which is not less than 0.2 megapascal at its normal operating temperature and which is currently not less than 0.2 megapascal, or compressed acetylene gas, the pressure of which is not less than 0.2 megapascal at a temperature of 15 degrees Celsius;

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(iv) In addition to what is listed in the preceding item, those liquefied gases, the pressure of which exceeds zero Pascal at a temperature of 35 degrees Celsius, and which, inclusive of liquefied hydrogen cyanide and liquefied methyl-bromide, are specified by a Cabinet Order.

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The gauge pressure is based on the atmospheric pressure (=zero), usually used in an engineering field. In physical and chemical fields, absolute pressure, which is based on vacuum, is used.

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The reinspection of containers shall not apply to containers with an internal volume or less 1 deciliter and to containers to be used without tight sealing.

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Security Regulation for General High-Pressure Gas

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Standard for storage space of containers and containers filled with high pressure gas

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Containers filled with high pressure must be always held below 40 °C

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The containers must be measured preventions for mechanical shock by fall or drop, and the damage of valve, and must not be handled roughly.

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### Technical standards for consumption

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### Regulation on Safety of Containers

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##### Methods of Labeling (Article 10)

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Labeling of the name of high pressure gas that may be filled in the container and properties

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"Flammable" in the case of Flammable Gas, and "Toxic" in the case of Toxic Gas

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Ultra Low-temperature container and the container filling a high pressure gas shall have undergone and passed the container re-inspection by the Minister of METI every below period.

密閉型の窒素容器(=超低温容器) Ultra Low-temperature container	気密・断熱・附属品等の試験 Test of gas leakage, thermal insulation, accessories
製造後の経過年数 The number of years since production	20年未満 5年毎 Less than 20 years Every 5 years
	20年以上 2年毎 More than 20 years Every 2 years

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
Comparative inspection of pressure gauge should be done ever year.

通常のボンベ Normal container(bomb)	外観・耐圧・質量試験 inspection of outside appearance, pressure proof test, mass test
	5年毎 Every 5 years

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The date of the previous inspection is punched in the container.

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The pressure regulator must be installed in the bomb when using high pressure gas from the bomb.  
It is very dangerous to open the valve of the bomb without the pressure regulator.

- 圧力調整器の選定 Selection of the pressure regulator

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The regulator fitted to the gas should be used because there is a material prohibited using according to a sort of gases.

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The pressure range of the regulator should be selected to fit the purpose of use.

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The regulator should be installed tightly while be careful about the right and left screw.

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Handling of regulator should be done accurately, slowly and certainly.

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Do not face the pressure gauge of the regulator when you open the valve of bomb.

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The regulator should be removed after closing valve and releasing the inner gas of regulator.

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For nitrogen gas (normal gases): Right screw

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### Notes for transporting the containers of Liq. N<sub>2</sub> and He by elevator

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When liquid nitrogen or liquid helium is transported by elevator, the containers should travel alone in the elevator, without any people accompanying them inside the elevator. Also, the containers should bear a warning notice of the prescribed format to prevent anyone from getting in the elevator with them.

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
Note 1) The prescribed format for the notice can be downloaded from the web page of the Natural Science Center for Basic Research and Development's Low Temperature Instrumental Analysis Department at

<http://www.sci.hiroshima-u.ac.jp/kiki/>, or from the web page of the Safety and Health Committee at <http://home.hiroshima-u.ac.jp/anzen/>.

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